

Oasis Velocity can improve your throughput and image quality

BRAIN	Oasis				Velocity				Spatial Resolution	Time Reduction
	Sequence	FOV	Sl.Thickness	Matrix	Scan Time	FOV	Sl.Thickness	Matrix		
Sag T1 FLAIR	220	5	320 x 256	3:43	230	5	320 x 256	2:09	Equivalent	40%
Axial T2 FSE	220	5	512 x 320	3:07	220	5	512 x 320	1:19	Equivalent	
Axial T2 FLAIR	220	5	256 x 256	3:34	220	5	288 x 256	2:30	Higher	
Axial DWI	255	5	128 x 128	1:09	240	5	128 x 128	1:01	Equivalent	
Total Study Time	11:33				6:59					

CERVICAL SPINE	Oasis				Velocity				Spatial Resolution	Time Reduction
	Sequence	FOV	Sl.Thickness	Matrix	Scan Time	FOV	Sl.Thickness	Matrix		
Sag T2 FSE	240	3	384 x 288	3:23	240	3	384 x 288	1:57	Equivalent	42%
Sag T1 FLAIR	240	3	256 x 256	3:35	240	3	256 x 256	2:08	Equivalent	
Sag T2 FSE FatSat	240	3	320 x 224	4:39	240	3	320 x 256	2:07	Higher	
Axial T2 FSE	200	3	256 x 224	4:12	200	3	256 x 224	2:58	Equivalent	
Total Study Time	15:49				9:10					

LUMBAR SPINE	Oasis				Velocity				Spatial Resolution	Time Reduction
	Sequence	FOV	Sl.Thickness	Matrix	Scan Time	FOV	Sl.Thickness	Matrix		
Sag T2 FSE	280	4	320 x 288	3:21	280	4	320 x 320	1:58	Higher	41%
Sag STIR	280	4	320 x 192	3:49	280	4	320 x 224	2:01	Higher	
Sag T1 FSE	280	4	320 x 288	3:48	280	4	320 x 288	2:11	Equivalent	
Axial T2 FSE	200	4	256 x 192	3:58	200	4	288 x 192	2:32	Higher	
Cor T2 FSE	280	4	320 x 256	2:57	280	4	320 x 256	1:52	Equivalent	
Total Study Time	17:53				10:34					

SHOULDER	Oasis				Velocity				Spatial Resolution	Time Reduction
	Sequence	FOV	Sl.Thickness	Matrix	Scan Time	FOV	Sl.Thickness	Matrix		
Sag T2 FSE	150	3.5	288 x 256	4:21	150	3.5	288 x 256	2:00	Equivalent	34%
Axial PD FSE FatSat	160	4	288 x 224	3:33	160	4	288 x 224	2:41	Equivalent	
Cor T2 FSE FatSat	150	3.5	256 x 192	4:21	150	3.5	320 x 224	3:16	Higher	
Cor T1 FSE	150	3.5	288 x 256	3:36	150	3.5	320 x 256	2:22	Higher	
Axial T1 FSE FatSat	160	4	288 x 224	3:56	160	4	320 x 224	2:37	Higher	
Axial T2* GE (In Phase)	160	4	288 x 192	3:33	160	4	320 x 192	2:21	Higher	
Total Study Time	23:20				15:17					

KNEE	Oasis				Velocity				Spatial Resolution	Time Reduction
	Sequence	FOV	Sl.Thickness	Matrix	Scan Time	FOV	Sl.Thickness	Matrix		
Axial PD FSE FatSat	160	4	288 x 224	3:55	150	4	320 x 256	2:14	Higher	42%
Sag PD FSE	150	3	320 x 320	3:55	150	3	384 x 320	2:04	Higher	
Sag T2 FSE	150	3	320 x 320	2:56	150	3	384 x 384	2:03	Higher	
Cor PD FSE FatSat	150	3	256 x 256	3:47	150	3	320 x 256	2:28	Higher	
Cor T1 FSE	150	3	256 x 256	3:49	150	3	320 x 288	1:51	Higher	
Total Study Time	18:22				10:40					

ANKLE	Oasis				Velocity				Spatial Resolution	Time Reduction
	Sequence	FOV	Sl.Thickness	Matrix	Scan Time	FOV	Sl.Thickness	Matrix		
Axial T1 FSE	150	4	256 x 256	3:04	150	4	256 x 256	2:04	Equivalent	30%
Axial T2 FSE FatSat	150	4	288 x 192	4:03	150	4	288 x 192	2:47	Equivalent	
Sag T2 FSE FatSat	150	3	224 x 192	3:53	150	3	256 x 192	2:42	Higher	
Cor T2 FSE FatSat	150	4	288 x 192	3:46	150	4	288 x 192	2:37	Equivalent	
Axial T1 FSE FatSat	150	4	256 x 192	4:12	150	4	256 x 192	2:58	Equivalent	
Total Study Time	18:58				13:08					

HIPS	Oasis				Velocity				Spatial Resolution	Time Reduction
	Sequence	FOV	Sl.Thickness	Matrix	Scan Time	FOV	Sl.Thickness	Matrix		
Bilateral Cor T1 FSE	360	4	384 x 256	3:28	360	4	384 x 256	2:35	Equivalent	29%
Bilateral Cor STIR	360	4	320 x 224	3:46	360	4	320 x 256	3:19	Higher	
Uni Axial PD FSE FatSat	200	4.5	256 x 192	4:47	180	4	256 x 192	3:03	Higher	
Uni. Sag PD FSE	200	4.5	256 x 200	4:00	180	4	288 x 224	2:27	Higher	
Uni. Cor PD FSE	200	4	320 x 192	3:48	180	4	320 x 192	2:44	Higher	
Uni. Axial T1 RSSG FatSep	200	4.5	224 x 192	4:30	180	4	224 x 192	3:01	Higher	
Total Study Time	24:19				17:09					

Velocity's RF coil technology and IP-RAPID 2D accelerated imaging drive reduced scan times and increased spatial resolution flexibility